

## CYTOCHROME OXIDASE PROTOCOL

### PRINCIPLE:

Cytochrome oxidase (COX) is the collective name for Complex IV of the oxidative respiratory chain of enzymes located in mitochondria. This method is a modification of the "Nadi" reaction. The use of 3,3' diaminobenzidine (DAB) results in a brown insoluble compound at the site of cytochrome oxidase activity.

### SPECIMEN REQUIRED:

Snap frozen human striated muscle

### Controls:

Stain several different muscles simultaneously. There is normally always enzyme activity present in muscle. The degree and location is interpreted.

### METHOD

**Fixation:** None, use snap frozen tissue.

**Technique:** Cut 10 - 16 micron (12  $\mu$ m) sections in cryostat from snap frozen biopsy. Attach one or more sections to a No.1½, 22 mm square coverslip.

### Equipment:

Ceramic staining rack - Thomas Scientific #8542-E40  
Columbia staining dish - Thomas Scientific #8542-C12  
Columbia staining dish(jar) - Thomas Scientific #8542-E30  
Forceps      Latex gloves

### Reagents:

Catalase (Sigma C-10)  
Cytochrome C (Sigma C-2506)  
deionized water  
3, 3' Diaminobenzidine tetrahydrochloride (Sigma D-5637)  
**CARCINOGENIC**, store at -20°C  
Permout - Fisher SP15-100, **FLAMMABLE HEALTH HAZARD**  
Reagent alcohol, ACS histochemical Fisher A962-4, or HPLC A995  
**FLAMMABLE, TOXIC**, store at room temp. in flammable cabinet  
Sodium dibasic phosphate (Na<sub>2</sub>HPO<sub>4</sub>) anhydrous ACS  
(FW 141.96)-Sigma S9763, Fisher S374, or Mallinckrodt 7917  
Store at room temperature  
Sodium monobasic phosphate (Na<sub>2</sub>HPO<sub>4</sub>) monohydrate (FW 137.99)  
Store at room temperature  
Sucrose -Sigma S7903, Store at room temperature  
Xylenes - Fisher #HC700-1GAL, **FLAMMABLE**,  
Store at room temperature in flammable cabinet)

### Solutions:

- I. 0.2 M Phosphate Buffer, pH 7.6 (at room temperature)  
0.2 M sodium dibasic phosphate ( $\text{Na}_2\text{HPO}_4$ ) anhydrous (28.39gm/liter) ~95 ml  
0.2 M sodium monobasic phosphate ( $\text{Na}_2\text{HPO}_4$ ) Monohydrate (27.8 gm/liter deionized  $\text{H}_2\text{O}$ ) ~11 ml  
Add the monobasic to the dibasic while adjusting the final pH to ~7.6 at room temperature.
- II. Incubating Solution:  
7 mg DAB  
12 mg cytochrome C  
20  $\mu\text{g}$  catalase (approximately a few crystals)  
750 mg Sucrose  
5 ml 0.2 M Phosphate Buffer  
5 ml deionized  $\text{H}_2\text{O}$   
Dissolve well.
- III. Alcohol 50 % (v/v)  
reagent alcohol ~50 ml  
deionized water ~50 ml
- IV. Alcohol 70 % (v/v)  
reagent alcohol ~70 ml  
deionized water ~30 ml
- V. Alcohol 80 % (v/v)  
reagent alcohol ~80 ml  
deionized water ~20 ml
- VI. Alcohol 95 % (v/v)  
reagent alcohol ~95 ml  
deionized water ~ 5 ml

### Staining Procedure

1. Place coverslips in the incubating solution in a columbia staining dish (Thomas Scientific #8542-C12) for at least 1 hour but no more than 2 hours at  $37^\circ\text{C}$  or **overnight at room temperature**.
2. Wash with three exchanges of tap or deionized  $\text{H}_2\text{O}$ .
3. Dehydrate in ascending alcohols (50%,70%,80%,95% x 2, 100% x 2) in Columbia staining dish(jar) - Thomas Scientific #8542-E30.
4. Clear with at least 2 changes of xylene in Columbia staining dish - Thomas Scientific #8542-E30
5. Mount with PERMOUNT or another synthetic organic mounting medium.

Washington University School of Medicine  
Neuromuscular Lab  
CAP: 1923316  
CLIA: 26D0652044  
NY: PFI 3499

## **Results**

Sites of cytochrome oxidase activity are colored brown.

## **REFERENCES**

1. Seligman et al., 1968. In: *HISTOCHEMISTRY - THEORETICAL and APPLIED*, A.G. Everson Pearse, Vol. 2, 3rd Edition; Williams & Wilkins, Baltimore, 1972.

12/26/2015