

## **ESTERASE STAINING: ALPHA-NAPHTHYL ACETATE PROTOCOL**

### **PRINCIPLE:**

This protocol demonstrates the sites of non-specific esterases in tissue sections. This modification of the technique described by B. J. Davis demonstrates small, angular denervated muscle fibers and neuromuscular junctions

### **SPECIMEN REQUIRED**

Snap frozen human striated muscle

### **METHOD**

**Fixation:** None. Use snap frozen tissue

**Technique:** Cut 10 -16 micron (12 micron) sections in cryostat from snap frozen biopsy. Attach one or more sections to a No. 1 ½ 22 mm, square coverslip

### **Equipment:**

Ceramic staining rack - Thomas Scientific #8542-E40  
Columbia staining dish - Thomas Scientific #8542-C12  
Columbia staining dish(jar) - Thomas Scientific #8542-E30  
Forceps      Latex gloves  
Filter paper (Baxter #f2217-070, Grade 363 Qualitative)

### **Reagents:**

Acetone - Baxter #010-4 **FLAMMABLE**  
Alpha-naphthyl acetate - Sigma N8505 - Store desiccated at -20°C  
deionized water  
Hydrochloric acid, ACS - Fisher A144-500, **CORROSIVE**,  
Store at room temperature  
Basic Fuchsin- Santa Cruz 203731 (troubleshooting: RC;30/172),  
Store at room temperature  
Permout - Fisher SP15-100, **FLAMMABLE HEALTH HAZARD**  
Reagent alcohol, ACS - histochemical Fisher A962-4 or HPLC A995,  
**FLAMMABLE, TOXIC, TERATOGENIC**  
Store at room temperature in flammable cabinet  
Sodium nitrite certified crystalline - Fisher S347 or Sigma S2252,  
**STRONG OXIDIZER, COMBUSTIBLE**  
Sodium dibasic phosphate (Na<sub>2</sub>HPO<sub>4</sub>) anhydrous, ACS (FW 141.96)  
Store at room temperature  
Xylenes - Fisher #HC700-1GAL, **FLAMMABLE**,  
Store room temperature in flammable cabinet)

**Solutions:**

1. 0.2 M Sodium Phosphate  
Sodium phosphate dibasic, anhydrous ( $\text{Na}_2\text{HPO}_4$ ) 7.1 g  
deionized water → 250 ml  
Store at room temperature
  
2. 4% Basic Fuchsin HCl  
Dissolve Basic Fuchsin 0.5 g  
in deionized water 10 ml  
Heat gently on a hot plate (**DO NOT BOIL**)  
Add concentrated (12N) hydrochloric acid 2.5 ml  
Cool to room temperature  
Filter (Baxter #f2217-070, Grade 363 Qualitative)  
Store refrigerated (0-5°C)
  
3. 4% Sodium Nitrite  
Sodium nitrite ( $\text{NaNO}_2$ ) 0.5 g  
deionized water 12.5ml  
Store refrigerated (0-5°C)
  
4. "Azotized Basic Fuchsin" (**PREPARED FRESH FOR EACH STAIN**)  
4% Basic Fuchsin -HCl (Solution #2) 0.2 ml  
4% sodium nitrite (Solution #3) 0.2 ml  
Sit at room temperature for a few minutes. (Amber color solution.)
  
5. Staining Solution (**PREPARED FRESH FOR EACH STAIN**)  
  
Into a 30 ml glass beaker **ADD IN THE ORDER STATED**  
Alpha-naphthyl Acetate ~ 2 mg  
Acetone ~ 0.75 ml  
MIX WELL  
Add 0.2 M sodium phosphate ( $\text{Na}_2\text{HPO}_4$ ) 12.5 ml  
(*SOLUTION MAY BECOME CLOUDY - OK!*)  
MIX WELL  
  
Add Solution #4 ("Azotized Basic Fuchsin")  
and *MIX WELL*  
**SOLUTION WILL CHANGE COLOR FROM YELLOW TO RED-  
BROWN IN LESS THAN FIVE (5) MINUTES !**
  
6. Alcohol 50 %  
reagent alcohol ~50 ml  
deionized water ~50 ml

7. Alcohol 70 %	
reagent alcohol	~70 ml
deionized water	~30 ml
8. Alcohol 80 %	
reagent alcohol	~80 ml
deionized water	~20 ml
9. Alcohol 95 %	
reagent alcohol	~95 ml
deionized water	~ 5 ml

### Staining Procedure

1. Place coverslips into a Columbia staining dish (Thomas Scientific #8542-012)
2. When the Staining Solution (Solution %) is gold to red-orange in color, add it to the coverslips in the staining dish for 5 minutes at room temperature.
3. Immediately place sections under running tap water for several minutes to wash the reaction product off the sections.
4. Clean back of coverslips with cotton swab.
5. Place coverslips with sections in a ceramic rack (Thomas Scientific #8542-E40).
6. Dehydrate in ascending alcohol solutions (50%, 70%, 80%, 95% x 2, 100% x2) - in Columbia staining dishes - Thomas Scientific #8542-E30.
7. Clear with xylene (x - 4x) also in Columbia staining dishes - Thomas Scientific #8542-E30.
8. Mount coverslip onto a labeled glass slide with Permount or other suitable organic mounting medium.

### Results

Esterase activity is demonstrated in denervated fibers as a red-brown color. Normal fibers exhibit a pale yellow to brown color. Neuromuscular junctions are demonstrated by a dark red-brown deposit on the edge of muscle fibers.

### REFERENCES

1. Thompson, **SELECTED HISTOCHEMICAL AND HISTOPATHOLOGICAL METHODS**, Charles C. Thomas, Springfield, IL, 1966.